



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## VACUUBRAND suction system coated with 2 Ltr. glass bottle

**€995.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The suction system is designed for connection to an existing external vacuum source such as a diaphragm pump or vacuum network and does not require an electrical supply.

The chemical-resistant, mechanical vacuum regulator ensures an optimum working vacuum. The regulator opens only when needed, limiting the gas load for the vacuum source.

Evaporation of the collected liquid is minimized. The guide tube integrated in the bottle head serves to reduce foam and aerosol formation and ensures long service life of the hydrophobic protective filter.

**Order No. 7-8155:** 2 liter coated glass bottle

### Performance Features

Mechanical vacuum regulator to minimize gas load on the vacuum network and provide optimum working vacuum. High stability, no electrical connection necessary for flexible use. Optimally equipped autoclavable collection bottle to reduce aerosol and foam formation, autoclavable quick disconnect couplings available as accessories. Suction handle with continuous suction hose - mechanics not in contact with media - no contamination. Hydrophobic protective filter 0.2 µm - high safety against contamination of the vacuum network, autoclavable as a unit with the collection bottle.

Weight: approx. 2.8 kg

Accessories: Spare collection bottle (Order No. 7-8156 and 7-8157)

Quick coupling complete (Order No. 7-8158)

## Additional Information

No.	VZ-1460
Manufacturer (Brand)	VACUUBRAND
old neoLab article no.	7-8155
EAN	4058072426750
Transport temperature	Room temperature
Wide	180 mm
Height	430 mm
Length	230 mm
Depth	230 mm
Weight	2.8 kg
Pumps and suction systems Features	with shatterproof coating with safety filter with collecting bottle
Functionality	mechanical
for DM hose inside 1	6.0 to 10.0 mm
Spare part	VZ-0115
Material collection bottle	Borosilicate glass
Type suction device/system	without pump
Type connection collection bottle	Hose olive
Volume collection bottle	2 l

