



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLabLine neoVAQ aspirator with vacuum pump, bottle 0.5 Ltr, LxHxW: 27x13x22 cm

€745.20
plus VAT &
Shipping

Product Images



Description

The neoVAQ Aspirator Portable is characterized by its compact design, which makes it easy to transport. It is ideally suited for work in cell culture, microbiology or when aspirating liquid waste in chemistry.

With matching slot for the aspirator.

Safety filter: 0.45 µm pore size (PTFE)

Complete with handle, 1-channel and 8-channel tip attachment for pipette tips with ejector. With stainless steel holder for 50 ml centrifuge tubes, with matching slot for the aspirator.

Product video:

Additional Information

No.	2-9345
Manufacturer (Brand)	neoLabLine
EAN	4058072085513
Transport temperature	Room temperature
Wide	265 mm
Height	210 mm
Length	130 mm
Depth	130 mm
Weight	1.3 kg
Voltage	240 V
Pumps and suction systems Features	with vacuum pump with safety filter with collecting bottle
Application pump	Evacuate / Suction
Pressure MIN	670.0000 mbar abs.
Spare part	2-9346 2-9347 1-0197
Functionality	electric
MAX flow rate [l/min]	1 l/min
Material collection bottle	Borosilicate glass
Type suction device/system	with vacuum pump
Type connection collection bottle	Hose olive
Type pump	Extraction devices and systems
Volume collection bottle	0.5 l
for DM hose inside l	5
for medium	Liquids Cell cultures Aqueous solutions Viscous media

