



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLabLine neoVAQ suction device with vacuum pump, bottle 4 Ltr

€1,144.00
plus VAT &
Shipping

Product Images



Description

The **neoVAQ aspirator Portable 4.0** combines compact design with high capacity. On the one hand, it is characterized by its compact design, which makes it easy to transport with the help of the carrying handle. At the same time, with 4 liters, it offers the usual capacity in terms of waste volume. The vacuum level and thus the suction force on the suction handle can be finely adjusted via a rotary knob. This ensures sensitive suction of liquid supernatants. The hoses can be easily connected to the lid of the 4-liter waste bottle made of polypropylene (PP) via quick-release couplings. In the lid of the 4 liter waste bottle there is a mechanical overflow protection, which prevents the extracted liquid from getting out of the bottle. In addition, the vacuum pump is protected by a safety filter made of PTFE. A plastic-coated glass bottle with olive connections is also available as an accessory.

It is ideal for work in cell culture, microbiology or when aspirating liquid waste in chemistry.

- 4 liter waste bottle made of PP, autoclavable
- Safety filter: 0.45 µm pore size (PTFE)
- with quick couplings
- with overflow protection in the bottle

Complete with hand operator, 1-channel and 8-channel tip attachment for pipette tips with ejector.

Additional Information

No.	2-9360
Manufacturer (Brand)	neoLabLine
EAN	4058072029067
Transport temperature	Room temperature
Wide	335 mm
Height	300 mm
Length	245 mm
Depth	245 mm
Weight	7 kg
Voltage	240 V
Pumps and suction systems Features	with vacuum pump with safety filter with collecting bottle
Application pump	Evacuate / Suction
Pressure MIN	30.0000 mbar abs.
Spare part	2-9361 2-9362 6-0035 6-0036
Functionality	electric
MAX flow rate [l/min]	20 l/min
Material collection bottle	Polypropylene (PP)
Type suction device/system	with vacuum pump
Type connection collection bottle	Hose olive Quick coupling
Type pump	Extraction devices and systems
Volume collection bottle	4 l
for DM hose inside l	5
for medium	Liquids Cell cultures Aqueous solutions Viscous media

