



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Bürkle dosing pump Dosi-Pump 250 ml,
immersion depth 41.5 cm, location hole Ø 28.5
mm, fixed outlet**

€107.90 As low
as
€81.60
plus VAT &
Shipping

Product Images

Description

Standard dosing approx. 250 ml per stroke, the pump has a fixed outlet. Dosing pumps suitable for drums, canisters, balloons or bottles, for dosing liquids of different viscosity. The dosing pumps are mounted by you in your existing container closure (see table mounting hole). For this purpose, an appropriate hole is drilled in the closure, the pump is inserted through it and screwed against it from below. A suitable container connection with location hole is not included in the scope of delivery. Wetted parts made of PP, PE and stainless steel

Additional Information

No.	KL-1102
Manufacturer (Brand)	Bürkle
VGKL number	100457250
Dual use indicator	false
EAN	4058072270230
Transport temperature	Room temperature
Net weight	0.17050 kg
sterile	No
Country of origin	Great Britain
Customs tariff number	84132000

