



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF Dosing pump SIMDOS 02 FEM1.02FT.18RCP

€3,443.00
plus VAT &
Shipping

Product Images

Description

Chemically resistant diaphragm metering pump SIMDOS® 02 FT.18 RC P - Delivery rate from 0.03 - 20 ml/min / pressure head 60 mWS, suction head 2 mWS - Pump head in PVDF material design - FFKM diaphragm, FFKM valves - Delivery rate adjustable manually and by means of external control as well as including RS 232 interface (RC Plus version) - Additional safety diaphragm for maximum protection - Easy replacement of the delivery diaphragm by activating the maintenance command stored in the operating program Technical data Application Dosing of liquids Delivery rate (ml/min) with water of 20 °C at zero delivery head 0,03 - 20 Pressure head (mWS) 60 Suction head (mWS) 2 Hose connection (mm) ID1.6 / AD3.2 Permissible media and ambient temperature Ambient temperature +5 ... +40 °C Media temperature +5 ... +80 °C Weight (kg) 0.9 Dimensions W x H x D (mm) 93 x 144 x150 Material pump head PTFE Diaphragm FFKM Valves FFKM Included in delivery: -Dosing pump -Suction filter -2 meter connection hose - Power supply unit Accessories (not included in scope of delivery) Stand bracket Order no. 160474 Wall bracket Order no. 160473 Foot switch for RCP version (RCP = manual operation and external control) Order no. 155872 In-line filter FS 60 T PVDF mesh size 70 µm Order no. 165210 FS 60 X PEEK mesh size 35 µm Order no. 165212 Pressure retaining valve JFDV30FTZ: 166288 (PTFE) Pressure retaining valve JFDV30TTZ: 166287 (PVDF) Pressure retaining valve JFDV30KTZ: 166286 (PP)

Additional Information

No.	KF-0436
Manufacturer (Brand)	KNF
VGKL number	106277467
EAN	4058072868123
Transport temperature	Room temperature

