



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle solvent pump Mini, 1.8 l/min, connection GL45

**€524.05**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

High quality solvent pump Mini for all common organic solvents such as acetone, isopropanol, ethanol, methanol etc.. Easy filling of small quantities from bottles and jugs up to 5 liters with GL 45 thread. Hose and seal made of chemically resistant PTFE.

Adapters to S 40, GL 38 and GL 32 available as accessories.

## Additional Information

---

No.	KL-1092
Manufacturer (Brand)	Bürkle
old neoLab article no.	2-2350
EAN	4058072270131
Transport temperature	Room temperature
Pumps and suction systems Features	with special chemical resistance
Application pump	Convey
Functionality	mechanical
Thread	GL 45
Material seals (wetted)	Polytetrafluoroethylene (PTFE)
Material pump (wetted)	Stainless steel
Type Connection Input	Screw connection
Type pump	Piston pumps. manual (barrel pump)
for barrels up to	5 l
for medium	Liquids

