



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF LIQUIPORT® NF 1.300 KT.18 RC diaphragm liquid pump

€1,136.00
plus VAT &
Shipping

Product Images

Description

Diaphragm liquid pump LIQUIPORT® NF 1.300 KT.18 RC - Delivery rate from 0.5 - 3 l/min / Pressure head 60 mWS, suction head 3 mWS - Self-priming, dry-running safe - Pump head in material version PP - Diaphragm made of PTFE, Valves made of FFKM - Flow rate adjustable manually and by means of external control (RC version) Technical data Application Pumping liquids Flow rate (l/min) for water at 20 °C at zero head 0.5 - 3 Pressure head (mWS) 60 Suction head (mWS) 3 Hose connection (mm) ID12 Permissible media and ambient temperature Ambient temperature +5 ... +40 °C Media temperature +5 ... +80 °C Weight (kg) 1.5 Dimensions W x H x D (mm) 104 x 188 x 160 Material pump head PP Diaphragm PTFE Valves FFKM Accessories Tripod bracket Order No. 160474 Wall bracket Order No. 160473 Foot slider for RC version (RC = manual operation and external control) Order No. 155872 Accessories (optional): Hose (material PTFE) 007860

Additional Information

No.	KF-0283
Manufacturer (Brand)	KNF
VGKL number	106273693
EAN	4058072868512
Transport temperature	Room temperature

