



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF LIQUIPORT® NF300KT.18S diaphragm liquid pump

€886.00
plus VAT &
Shipping

Product Images

Description

Diaphragm liquid pump LIQUIPORT® NF 300 KT.18 S - Delivery rate from 0.5 - 3 l/min / Pressure head 10 mWS, suction head 3 mWS - Self-priming, dry-running safe - Pump head in material version PP - Diaphragm made of PTFE, Valves made of FFKM - Pressure head also available for 60 mWS on request - Delivery rate manually adjustable (version S) Technical data Application Delivery of liquids Delivery rate (l/min) for water at 20 °C at zero head 0.5 - 3 Pressure head (mWS) 10 Suction head (mWS) 3 Hose connection (mm) ID12 Permissible media and ambient temperature Ambient temperature +5 ... +40 °C Media temperature +5 ... +80 °C Weight (kg) 1.5 Dimensions W x H x D (mm) 104 x 188 x 160 Material pump head PP Diaphragm PTFE Valves FFKM Accessories Tripod bracket order no. 160474 Wall bracket order no. 160473 Accessories (optional): Hose (material PE) 007859 Hose (material PTFE) 007860 Mounting set (stand bracket) 160474 Mounting set (wall bracket) 160473

Additional Information

| | |
|-----------------------|------------------|
| No. | KF-0274 |
| Manufacturer (Brand) | KNF |
| VGKL number | 566260020 |
| EAN | 4058072868420 |
| Transport temperature | Room temperature |

