



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND rotary vane pump RZ 2.5 + FO + V16 (package), completely assembled, ready for operation

Price on
request
plus VAT &
Shipping

Product Images



Description

The two-stage RZ 2.5 is a powerful rotary vane pump with a compact design, low weight and quiet running noise. It is ideal for laboratory and plant applications where a very good ultimate vacuum must be achieved at medium gas flow rates (chemistry, life sciences and physics).

Scope of delivery:

Pump completely assembled, ready for operation, with oil filling and instructions. With exhaust filter FO and valve VS 16.

Additional Information

No.	VZ-1361
Manufacturer (Brand)	VACUUBRAND
old neoLab article no.	7-0213
EAN	4066292393936
Transport temperature	Room temperature
Wide	210 mm
Height	190 mm
Length	316 mm
Depth	316 mm
Weight	12.2 kg
Volume	43 dB
Voltage	230 V
Pumps and suction systems Features	with special chemical resistance with valve with gas ballast with exhaust filter
Application pump	Evacuate / Suction
Pressure MIN	0.0020 mbar abs.
Functionality	electric mechanical
MAX flow rate [l/min]	31 l/min
Power	180 W
Number of steps	2
Protection class	IP 40
Type Connection Output	Hose olive
Type Connection Input	Hose olive
for DM hose inside l	10
Type pump	Rotary vane pumps
for medium	Air/Gases/Vapors
Oil filling MAX	280 ml
Oil filling MIN	120 ml

