



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® Mini rotary vane pump 4l/min.

**€266.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Small, low-cost **rotary vane pump for non-toxic, non-flammable and non-explosive gases**, operating completely oil-free according to the rotary principle.

The vanes are made of special carbon with service lives of up to 10,000 operating hours (in the partial load range).

The pump housing is made of stainless steel 1.4305.

- oil-free gases, no valves
- largely corrosion resistant to aggressive gases
- suitable for use in extreme climatic conditions from -30°C to +65°C
- Driven by powerful DC motors, with stranded wires, ball bearings, carbon brushes, 2 hose olives 4 mm (nickel-plated brass)
- high performance, low weight

## Additional Information

No.	1-1500
Manufacturer (Brand)	neoLab
EAN	4058072071257
Wide	30 mm
Height	25 mm
Length	59.6 mm
Depth	59.6 mm
Weight	98 g
Volume	70 dB
Voltage	6 V
Pumps and suction systems Features	oil-free
Application pump	Evacuate / Suction
Number of steps	1
Pressure MAX	160.00 bar abs.
Pressure MIN	50.0000 mbar abs.
Functionality	electric
MAX flow rate [l/min]	4 l/min
Power	5 W
Type Connection Output	Hose olive
Type Connection Input	Hose olive
Type pump	Mini rotary vane pump
for DM hose inside 1	4
for medium	Air/Gases/Vapors

