



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wiggens rotary vane pump, 50 L/min

**€833.95**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

oil-sealed **rotary vane vacuum pumps** have become indispensable in research and production.

They can be used independently and can be used as backing pumps for molecular pumps or diffusion pumps. Wiggins' direct-drive rotary vane pumps are available in single- and two-stage versions and cover the spectrum from desktop laboratory pumps to pumps for the production line.

Low noise and stable operating temperature for common voltages and frequencies, long life and excellent reliability for chemical synthesis.

Various applications can include vacuum distillation, vacuum filtration, vacuum testing, vacuum freeze drying, vacuum coating, vacuum leak detection, vacuum packaging, vacuum adsorption, etc.

## Additional Information

No.	WI-1517
Manufacturer (Brand)	Wiggins
old neoLab article no.	F-1000
EAN	4058072163266
HSNumber	8414100090
Transport temperature	Room temperature
Weight	7 kg
Volume	59 dB
Voltage	230 V
Pumps and suction systems Features	with special chemical resistance with gas ballast
Application pump	Evacuate / Suction
Pressure MIN	0.0250 mbar abs.
Functionality	electronic
MAX flow rate [l/min]	50 l/min
Power	180 W
Number of steps	1
Type Connection Output	Screw connection
Type Connection Input	Screw connection
Type pump	Rotary vane pumps
for medium	Air/Gases/Vapors
Oil filling MAX	250 ml

