




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle PP barrel pump with discharge bend, 80 cm, 260 ml/ stroke

€165.15 As low
as
€142.80
plus VAT &
Shipping

Product Images

Description

Drum pump PP with rigid discharge elbow. Barrel pumps made of polypropylene enable fast filling of liquids from barrels, canisters and tanks. Due to the good chemical resistance of PP, acids, alkalis and aggressive cleaners can also be filled. The deep bottom check valve of the barrel pump and the continuously adjustable immersion depth by means of the barrel screw connection ensure good residual emptying. The barrel screw connection secures the barrel pump in the barrel. Barrel screw connections and thread adapters for commercially available containers are available as accessories. Seals made of PTFE Immersion tube Ø 32 mm Delivery rate up to approx. 12 l/min For drums and tanks from 30-1000 liters Attention! "List of chemical resistances" must be observed!

Additional Information

No.	KL-1064
Dual use indicator	false
Manufacturer (Brand)	Bürkle
VGKL number	290510801
EAN	4058072269852
Transport temperature	Room temperature
Net weight	1.33950 kg
Country of origin	Germany
sterile	No
Customs tariff number	84132000

