



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle lever pump, yellow, stainless steel piston rod

**€59.35**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The lever pump is made of polypropylene with a piston rod made of stainless steel V2A (1.4301) and seals made of FKM. This manual lever pump can be used to remove quantities as required from drums without much effort. Ideal for occasional use. Applications of the lever pump: water, lubricants, light oils, as well as water-based chemicals and light acids and alkaline solutions. With integrated barrel screw connection for 1 1/4" and 2" fine thread Delivery rate 350 ml per stroke Immersion depth infinitely variable from 65-90 cm "Note list of chemical resistances"!

## Additional Information

---

No.	KL-1110
Dual use indicator	false
Manufacturer (Brand)	Bürkle
EAN	4058072270315
Transport temperature	Room temperature
Net weight	0.97000 kg
Country of origin	Taiwan, Province of China
sterile	No
Customs tariff number	84132000

