



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle magnetic centrifugal pump, 15 W, with hose nozzle

**€514.45** As low  
as  
**€437.55**  
plus VAT &  
Shipping

### Product Images



## Description

---

Magnetic centrifugal pump 15 Watt, motor 230 V, materials PP (glass fiber reinforced), stainless steel (shaft), Rulon (bearing), FKM (seal). The magnetic centrifugal pumps are suitable for pumping acids, alkalis, alkaline and galvanic baths, photographic chemicals as well as radioactive, sterile and highly corrosive liquids. Due to the drive via non-contacting magnetic force transmission, no dynamic shaft seals are required, thus the pump housing is hermetically sealed and does not allow any leakage. Flow rate 1 to 12 l/min, delivery head up to 2 m Without thread, hose nozzle 14 mm Sealless Ambient temperature max. 45 °C Operating temperature max. 60 °C Viscosity max. 20 mPas

## Additional Information

---

No.	KL-1140
Dual use indicator	false
Manufacturer (Brand)	Bürkle
EAN	4058072270612
Transport temperature	Room temperature
Net weight	1.13950 kg
Country of origin	Italy
sterile	No
Customs tariff number	84137051

