



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Manual barrel pump made of PP

€113.80
plus VAT &
Shipping

Product Images



Description

Manual drum pump Suitable for chemicals, food, water, AdBlue®, as well as for high viscosity liquids such as lubricants and cooking oils. Completely made of polypropylene, food safe. The capacity per stroke is approx. 400 ml.

Additional Information

No.	2-2309
Manufacturer (Brand)	neoLab
EAN	4058072060480
Transport temperature	Room temperature
Length	100 cm
Pumps and suction systems Features	with special chemical resistance
Application pump	Conveying dosing/transferring
Functionality	mechanical
Conveying capacity	0.4 l/stroke
Thread	BSP 2 inch
Material pump (wetted)	Polypropylene (PP)
Type Connection Input	Screw adapter
Material seals (wetted)	Fluoropolymer rubber (FPM/FKM)
Type pump	Piston pumps. manual (barrel pump)
for barrels up to	229 l
for medium	Organic solutions Aqueous solutions Oils Viscous media

