



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Syphon pump w. screw adapter, suction tube 90 cm

€46.00
plus VAT &
Shipping

Product Images



Description

Ideal for removing and transferring many oils and liquids.

The medium comes into contact only with polypropylene, stainless steel and HDPE. Simply immerse the suction tube and pump a few times with the suction piston. The syphonic effect can be enhanced by the suction piston.

No clogging, no clogging of the valve system.

Additional Information

No.	2-2303
Manufacturer (Brand)	neolab
EAN	4058072060442
Transport temperature	Room temperature
Length	90 cm
Pumps and suction systems Features	with special chemical resistance
Application pump	Conveying dosing/transferring
Functionality	mechanical
MAX flow rate [l/min]	19 l/min
Thread	Fine thread
Material pump (wetted)	Polypropylene (PP) Polyethylene (PE)
Type Connection Input	Screw connection
Material seals (wetted)	High-density polyethylene (HDPE)
Type pump	Syphon pumps
for medium	Oils Liquids

