



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle barrel pump PP, 300ml/stroke, R2" u. 3/4" outside, 84cm

**€36.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The piston pump for aqueous and petrochemical products. It offers a good price-performance ratio. The piston pump is mounted on the filling container by means of a 2" or 3/4" barrel thread. After mounting, the pump can still be rotated around its axis, so containers can also be placed on the drum lid for filling. Polypropylene, polyethylene, seals made of NBR Flow rate adjustable to 100, 250, 350 and 500 ml per stroke For containers with G 2" and G 3/4" internal threads Immersion tube length 84 cm "Note the list of chemical resistances"!

## Additional Information

---

No.	KL-1108
Manufacturer (Brand)	Bürkle
Dual use indicator	false
EAN	4058072270292
Transport temperature	Room temperature
Net weight	0.38000 kg
sterile	No
Country of origin	Taiwan, Province of China
Customs tariff number	84132000

