




neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wiggens oil-free piston pump, vacuum 100 mbar, 80 L/min

€962.24
plus VAT &
Shipping

Product Images



Description

Wiggens**oil-free piston vacuum pumps** with innovative electronic and mechanical technology. With modern design concept, they feature the compact design and light weight. Maintenance-free, safe and convenient.

Features

No air pollution.

Vacuum pumps operate without lubricants, regular oil changes and maintenance. No oil pollution.

Moisture trap with filter cartridge

The piston pumps are equipped with a filter cartridge in the air inlet to filter particles and moisture, and to extend the life of the pump.

Vacuum regulator

The oil-free piston pumps are equipped with a vacuum regulator to adjust the vacuum.

Thermal protection device

Each motor is equipped with a built-in thermal protection device to automatically shut down the pump in case of overheating and turn it on again after the temperature cools down.

Applications

- Biological laboratories
- Food industry
- Microbiological detection
- Vacuum extraction
- Liquid filtration
- Vacuum drying
- Measurement of suspended solids

Additional Information

No.	WI-1536
Manufacturer (Brand)	Wiggins
old neoLab article no.	F-1019
EAN	4058072163457
Transport temperature	Room temperature
autoclavable at 121°C	No
Wide	340 mm
Height	195 mm
Depth	185 mm
Weight	8.6 kg
Volume	52 dB
Voltage	230 V
Pumps and suction systems Features	oil-free
Application pump	Evacuate / Suction
Pressure MIN	100.0000 mbar abs.
Functionality	electric
MAX flow rate [l/min]	80 l/min
Power	220 W
Number of steps	1
Type pump	Piston pumps
for medium	Air/Gases/Vapors

