



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wiggens oil-free piston pump, vacuum 30 mbar, 50 L/min

**€962.24**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Wiggens**oil-free piston vacuum pumps** with innovative electronic and mechanical technology. With modern design concept, they feature the compact design and light weight. Maintenance-free, safe and convenient.

### Features

No air pollution.

Vacuum pumps operate without lubricants, regular oil changes and maintenance. No oil pollution.

### Moisture trap with filter cartridge

The piston pumps are equipped with a filter cartridge in the air inlet to filter particles and moisture, and to extend the life of the pump.

### Vacuum regulator

The oil-free piston pumps are equipped with a vacuum regulator to adjust the vacuum.

### Thermal protection device

Each motor is equipped with a built-in thermal protection device to automatically shut down the pump in case of overheating and turn it on again after the temperature cools down.

### Applications

- Biological laboratories
- Food industry
- Microbiological detection
- Vacuum extraction
- Liquid filtration
- Vacuum drying
- Measurement of suspended solids

## Additional Information

No.	WI-1537
Manufacturer (Brand)	Wiggins
old neoLab article no.	F-1020
EAN	4058072163464
Transport temperature	Room temperature
autoclavable at 121°C	No
Wide	340 mm
Height	195 mm
Depth	185 mm
Weight	8.6 kg
Volume	52 dB
Voltage	230 V
Pumps and suction systems Features	oil-free
Application pump	Evacuate / Suction
Pressure MIN	30.0000 mbar abs.
Functionality	electric
MAX flow rate [l/min]	50 l/min
Power	220 W
Number of steps	1
Type pump	Piston pumps
for medium	Air/Gases/Vapors

