



KNF Compressor N022AN.18 IP44

Product Images

Description

Diaphragm compressor pump N022AN.18 only for positive pressure applications (e.g. positive pressure filtration)
===== For vacuum applications (e.g. vacuum filtration, suction tasks, degassing, ...) we recommend the LABOPORT N816.3KT.18 pump -higher flow rate -significantly better vacuum -chemical-resistant, corrosion-resistant -narrow, space-saving design -quiet operation

Additional Information

No.	KF-0381
Manufacturer (Brand)	KNF
EAN	4058072866822

