



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND diaphragm pump MD 1, three-stage, 230 V / 50-60 Hz, CEE power cable

€1,700.00
plus VAT &
Shipping

Product Images



Description

Diaphragm pump MD 1 --- --- Max. Pumping speed 50/60 Hz (m³/h): 1.2/1.4 --- Ultimate vacuum (abs.) (mbar): 1.5 --- Connection suction side: hose shaft DN 6 mm --- Connection discharge side: silencer / G1/8" --- Dimensions (L x W x H), approx. (mm): 303x143x163 --- Weight, approx. (kg): 6.5 --- Nominal voltage / mains frequency: 200-230 V/50-60 Hz --- Plug: CEE --- Motor nominal power (kW): 0.08 --- Scope of delivery: pump completely assembled, ready for operation, with instructions. Old item number: 696080

Additional Information

No.	VZ-1256
Manufacturer (Brand)	VACUUBRAND
VGKL number	101871069
EAN	4066292392922
Gross weight	8.7
Transport temperature	Room temperature
Net weight	6.70000 kg
Country of origin	Germany
Customs tariff number	84141089

