



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND MD 4 NT diaphragm pump, three-stage, certification (NRTL): C/US 230 V / 50-60 Hz

€2,925.00
plus VAT &
Shipping

Product Images



Description

Diaphragm pump MD 4 NT --- --- Max. Pumping speed 50/60 Hz (m³/h): 3.8/4.3 --- Ultimate vacuum (abs.) (mbar): 1 --- Connection suction side: Small flange KF DN 16 --- Connection discharge side: Silencer / G1/4" --- Dimensions (L x W x H), approx. (mm): 325x239x198 --- Weight, approx. (kg): 16.4 --- Nominal voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Motor nominal power (kW): 0.25 --- Scope of delivery: Pump completely assembled, ready for operation, with instructions. Old item number: 736000

Additional Information

No.	VZ-1598
Manufacturer (Brand)	VACUUBRAND
EAN	4066292396227
Gross weight	19.7
Transport temperature	Room temperature
Net weight	17.20000 kg
Country of origin	Germany
Customs tariff number	84141089

