



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WELCH diaphragm pump MP 301 E 230V 50/60Hz

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Suction capacity 50/60Hz 38/41 l/min Motor power 180 W Suction connection hose shaft DN 8 for hose inner diameter 8 mm Pressure connection silencer and hose shaft DN 8 for hose inner diameter 8 mm (as accessory) Option gas ballast valve

Additional Information

No.	WL-0214
Manufacturer (Brand)	WELCH
VGKL number	109241171
EAN	4066292439696
Transport temperature	Room temperature

