



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND diaphragm pump MZ 2 NT, two-stage, 230 V / 50-60 Hz, CEE power cable certification

€1,540.00
plus VAT &
Shipping

Product Images



Description

Diaphragm pump MZ 2 NT --- --- Max. Pumping speed 50/60 Hz (m³/h): 2.2/2.4 --- Ultimate vacuum (abs.) (mbar): 7 --- Connection suction side: hose shaft DN 10 mm --- Connection discharge side: silencer / G1/4" --- Dimensions (L x W x H), approx. (mm): 243x239x198 --- Weight, approx. (kg): 11.0 --- Nominal voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Motor nominal power (kW): 0.18 --- Scope of delivery: pump completely assembled, ready for operation, with instructions. Old item number: 732000

Additional Information

No.	VZ-1497
Manufacturer (Brand)	VACUUBRAND
VGKL number	101871073
EAN	4066292395220
Gross weight	13.7
Transport temperature	Room temperature
Net weight	11.30000 kg
Country of origin	Germany
Customs tariff number	84141089

