




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF Mini Diaphragm Vacuum Pump N816.1.2KT.18

€881.00
plus VAT &
Shipping

Product Images



Description

Mini diaphragm vacuum pump N 816.1.2 KT.18 - Delivery rate 1.8 m³/h / Ultimate vacuum 160 mbar abs. - PTFE-coated diaphragm for aggressive/corrosive gases and vapors Technical data Application Filtration Delivery rate (m³/h) at atm. pressure 1.8 Ultimate vacuum (mbar abs.) 160 Pressure 1.8 Ultimate vacuum (mbar abs.) 160 Operating overpressure (bar) 0.5 Hose connections (mm) ID6 Permissible media and ambient temperature +5 ... +40 °C Weight (kg) 3.95 Dimensions W x H x D (mm) 102 x 141 x 361 Material Pump head PPS Diaphragm PTFE-coated Valves FFPM Accessories (optional): Fine adjustment valve with vacuum gauge (suction side): 057830 Noise damper: 000346 (mounting on pump head) Noise damper: 056366 (mounting on discharge side by means of hose) Hose nipple ID8: 004975 Vacuum chemistry hose ID 6mm (delivery in whole meters): 055535

Additional Information

No.	KF-0152
Manufacturer (Brand)	KNF
VGKL number	566210146
EAN	4058072867447
Gross weight	4.65
Transport temperature	Room temperature
Net weight	3.95000 kg
Country of origin	Germany
Customs tariff number	84148022
Voltage	230 V
Power	100 W
Protection class	IP 20

