



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KNF Vacuum Pump LABOPORT® N 96

€468.00
plus VAT &
Shipping

Product Images



Description

KNF Vacuum Pump, LABOPORT®, Type N96 (successor pump to N86KN.18, N86KT.18, N811KN.18, N811KT.18) The chemically resistant diaphragm pump Type N 96 is the smallest pump in our new LABOPORT® series and saves valuable space in the laboratory. The pump combines high flow rate with a compact design. The manual speed control makes the N 96 very energy efficient and quiet in use. The pumps of the new LABOPORT® generation are equipped with an easy-to-clean surface. Applications: filtration, SPE, liquid aspiration via vacuum, among others. Mini diaphragm vacuum pump LABOPORT®, type N 96: - Delivery rate 0.3 - 0.42 m³/h / ultimate vacuum < 130 mbar abs. - Very small footprint - Integrated speed control - PTFE-coated diaphragm for aggressive/corrosive gases and vapors - Included with the pump: 1x hose nozzle 6mm, 1x silencer/suction filter. 24V power supply with change-over plug Accessories Fine adjustment valve with vacuum gauge (vacuum): 132254/0000 Fine adjustment valve with pressure gauge (pressure): 058970 Hose nipple ID 6 mm 055956 Hose nipple ID 8mm 004975 Noise damper/suction filter 000346 Stand mounting 323484

Additional Information

No.	KF-0310
Manufacturer (Brand)	KNF
VGKL number	106270096
EAN	4058072868802
Transport temperature	Room temperature

