




Axygen® Money Documentation System

Product Images



Description

Money documentation system GD-1000 **for connection to a PC or notebook**.

The money documentation system provides 16 bit TIFF image files in publication quality. The commissioning of the device is simple and fast. The operation is easy and intuitive. Image files can be easily transferred to standard programs for gel evaluation.

Automatic exposure and capture. High resolution camera (5.4 MP). Safety switches prevent the UV lamp from shining when the door is open. The transillumination lamp has a lifetime of 30,000 hours, the user can choose between UV 302, UV 365 or Epi White.

Additional Information

No.	CQ-0380
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2032
EAN	4058072225797
Wide	381 mm
Length	304 mm
Operation via	externer PC
Voltage	230 V
Type interface	USB
Optical devices properties	mit integrierter Kamera
Energy supply	Stromnetz
Type illuminant	UV-Röhre
Wavelength	302 nm 470 nm
Light color	weiß
Model camera	integrierte 5.4 Megapixel Digitalkamera

