



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Axygen® Money Documentation System

**€6,719.29**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Money documentation system GD-1000 **for connection to a PC or notebook** .

The money documentation system provides 16 bit TIFF image files in publication quality. The commissioning of the device is simple and fast. The operation is easy and intuitive. Image files can be easily transferred to standard programs for gel evaluation.

Automatic exposure and capture .High resolution camera (5.4 MP). Safety switches prevent the UV lamp from shining when the door is open. The transillumination lamp a lifetime of 30,000 hours, the user can choose between UV 302, UV 365 or Epi White.

## Additional Information

---

No.	CQ-0380
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2032
EAN	4058072225797
Transport temperature	Room temperature
Wide	381 mm
Length	304 mm
Optical devices properties	with integrated camera
Operation via	external PC
Energy supply	Power grid
Light color	white
Model camera	integrated 5.4 megapixel digital camera
Voltage	230 V
Type illuminant	UV tube
Type interface	USB
Wavelength	302 nm 470 nm

