





Axygen® Money Documentation System BL

Product Images





Description

Standalone money documentation system GDBL-1000 including tablet .

The money documentation system provides 16 bit TIFF image files in publication quality. The commissioning of the device is simple and fast. The operation is easy and intuitive. Image files can be easily transferred to standard programs for gel evaluation.

Automatic exposure and capture. High resolution camera (5.4 MP). Safety switches prevent the UV lamp from shining when the door is open. The transillumination lamp a lifetime of 30,000 hours, the user can choose between UV 302, UV 365, Epi White and Epi Blue.

Additional Information

No.	CQ-0382
Manufacturer (Brand)	Axygen
old neoLab article no.	C-2033
EAN	4058072225810
Wide	381 mm
Length	304 mm
Operation via	Tablet Computer
Voltage	230 V
Type interface	USB Bluetooth
Optical devices properties	mit integrierter Kamera
Energy supply	Stromnetz
Type illuminant	UV-Röhre LED
Wavelength	302 nm 366 nm 470 nm
Light color	blau weiß
Model camera	integrierte 5.4 Megapixel Digitalkamera

