



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® LED light panel, light area: 42.9 x 30.9 cm

**€303.60
plus VAT &
Shipping**

Product Images



Description

Super flat, dimmable neoLab LED light plate with integrated rechargeable battery.

For uniform and color-correct illumination when viewing and sorting gels, plates, slides and other transparent specimens. Various power supply options: via the built-in lithium-ion battery, the supplied mains/charger or with the USB cable, e.g. via power bank, PC or car. Color temperature approx. 5000 Kelvin and color rendering index CRI=95 for correct color reproduction, both when judging and photographing originals. cm and inch scales next to the illuminated area.

A short touch of the touch switch activates the illuminated surface, a longer touch results in a stepless change of brightness. Scope of delivery: light panel, power supply/charger, USB cable

Additional Information

No.	2-2122
Manufacturer (Brand)	neoLab
EAN	4058072189327
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	black white
Material	Acrylonitrile butadiene styrene (ABS)
Wide	350 mm
Height	8 mm
Length	499 mm
Depth	499 mm
Weight	1.69 kg
Voltage	3 V
Energy supply	Battery
Power	7 W
Type illuminant	LED
Type battery	Battery

