



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Polarimeter OAB-N

€920.00
plus VAT &
Shipping

Product Images



Description

The ideal helper when starting to analyze your optically active solutions in the laboratory

- The Kern Polarimeter is a manual polarimeter that features ergonomic design and ease of use
- The powerful 589 nm sodium vapor lamp is an optimal light source for generating a linearly polarized light beam
- Precise determination of the angle of rotation of the substance under investigation is made possible by the 1° scale division including vernier (0.05°)
- Two glass cuvettes (100 mm/200 mm) are also included in the delivery for ideal recording of liquid samples
- Included in delivery: Sodium vapor lamp, 100 mm glass cuvette, 200 mm glass cuvette, spare lenses and sealing rings for cuvettes

Technical data

- Light source: Sodium vapor lamp (589 nm)
- Stabilization time: 10 min after switching on
- Overall dimensions W×D×H 500×130×330 mm

Main application areas

- Pharmacies
- Hospitals
- Beverage industry
- Food industry
- Chemical industry
- Laboratories
- Education

Additional Information

No.	KP-9859
Manufacturer (Brand)	KERN
EAN	4058072855208
Transport temperature	Room temperature
Color	White
Wide	500 mm
Height	330 mm
Number of lenses	1
Number of places objective turret	1
Number of scales	1
Operation via	Knobs/Buttons
Energy supply	Power grid
Weight	4.7 kg
Manufacturer	CORE
Voltage	230 V
Type lighting	Transmitted light
Type illuminant	Fluorescent lamp

