



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® handheld refractometer for antifreeze control and battery acid test

€213.00
plus VAT &
Shipping

Product Images



Description

Enables **measurement of the freezing point of ethylene and propylene glycol solutions.**

Particularly suitable for checking the state of charge of lead-acid batteries and for monitoring the concentration of heat exchange fluid in refrigeration systems.

Additional Information

No.	7-0114
Manufacturer (Brand)	neoLab
EAN	4058072027117
Transport temperature	Room temperature
Wide	40 mm
Height	29 mm
Length	190 mm
Depth	190 mm
Weight	0.22 kg
Optical devices properties	with automatic temperature compensation (ATC)
Number of scales	2
Refractometer measures	Alcohol content
Measured variables salinity	sg
Precision	5 °C 0.01 sg
Refractometer is suitable for	Battery fluids Heat exchange fluids
Scale ranges Salinity	0 to -50 °C 1.15 to 1.30 sg
Refractometer type	Hand refractometer

