



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Handheld refractometer for salt concentration, 0-28% Brix Salinity

€177.80
plus VAT &
Shipping

Product Images



Description

Refractometer **for determining the proportion of sodium chloride in seawater** and in solutions of prepared foods for a salt concentration of up to 28 %, i.e. 28 g salt in 100 g solution.

Also suitable for checking brines and for pickling liquids (canning industry).

The instrument is calibrated to the refractive index of sodium chloride.

Additional Information

No.	7-0110
Manufacturer (Brand)	neolab
EAN	4058072027094
Transport temperature	Room temperature
Wide	40 mm
Height	27 mm
Length	160 mm
Depth	160 mm
Weight	0.2 kg
Number of scales	1
Refractometer measures	Salinity
Measured variables salinity	% Brix Salinity
Precision	0.2% Brix Salinity
Refractometer is suitable for	Sea water Solutions of prepared food with salt content Salt brines
Scale ranges Salinity	0 to 28% Brix Salinity
Refractometer type	Hand refractometer

