



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Handheld refractometer for weak salt concentration, 0-100‰

**€177.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Measures salinity in weak concentration on a thousandths scale ‰

This easy-to-use refractometer quickly provides direct measurement of salinity and specific gravity of seawater.

## Additional Information

---

No.	7-0111
Manufacturer (Brand)	neoLab
EAN	4058072027100
Transport temperature	Room temperature
Wide	40 mm
Height	27 mm
Length	190 mm
Depth	190 mm
Weight	0.2 kg
Number of scales	1
Measured variables salinity	sg
Precision	0.001 sg
Refractometer is suitable for	Seawater
Refractometer measures	Salinity
Scale ranges Salinity	1.000 to 1,070 sg 0 to 100 ‰
Refractometer type	Hand refractometer

