



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN refractometer ORA,,468-1,518

€360.00
plus VAT &
Shipping

Product Images



Description

- The models of the KERN ORA-E series are universal, maintenance-free analog handheld refractometers
- The handy and robust design allows an easy, efficient and durable use in everyday life
- The manual conversion effort is avoided by several selectable scales and excludes application errors
- These scales are specially developed, precisely calculated and verified. They are also characterized by very thin and clear lines
- The optical system and the prism cover are made of special materials, which allow a low-tolerance measurement
- All models are equipped with an eyepiece that can be easily and smoothly adjusted to different visual acuities
- Models marked with "ATC" feature automatic temperature compensation, which enables accurate measurements at different ambient temperatures

Included in delivery:

- Storage box
- Calibration solution
- Pipette
- Screwdriver
- Cleaning cloth
- Further accessories are optionally available

Range of application Expert applications

The following models have a special large measuring range for refractive index and large split scales for measuring Brix values.

Main application areas:

- Universal range of use, especially in applications with a need for an extra large measuring range

Additional Information

| | |
|-----------------------|--------------------|
| No. | KP-5056 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072146092 |
| Transport temperature | Room temperature |
| Manufacturer | CORE |
| Refractometer type | Hand refractometer |

