



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN refractometer ORA,urine (sG) 1.000-1.061

€90.00
plus VAT &
Shipping

Product Images



Description

- The models of the KERN ORA-P series are universal, maintenance-free analog handheld refractometers
- The handy and robust design allows an easy, efficient and durable use in everyday life
- The manual conversion effort is avoided by several selectable scales and excludes application errors
- These scales are specially developed, precisely calculated and verified. They are also characterized by very thin and clear lines
- The optical system and the prism cover are made of special materials, which allow a low-tolerance measurement
- All models are equipped with an eyepiece that can be easily and smoothly adjusted to different visual acuities
- Models marked with "ATC" feature automatic temperature compensation, which enables accurate measurements at different ambient temperatures

Included in delivery:

- Storage box
- Calibration solution
- Pipette
- Screwdriver
- Cleaning cloth
- Further accessories are optionally available

Urine application range

The following models are particularly suitable for measuring specific urine weight (density), serum content (serum protein in urine) and refractive index.

Main application areas:

- Hospitals
- Medical practices
- Medical training facilities
- Homes for the elderly and nursing homes
- Sports medicine (doping control)
- Veterinarians

Additional Information

No.	KP-5083
Manufacturer (Brand)	KERN
EAN	4058072146368
Transport temperature	Room temperature
Manufacturer	CORE
Measured variables salinity	sg
Refractometer measures	Salinity
Scale ranges Salinity	1.000 to 1,060 sg
Refractometer type	Hand refractometer

