




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN refractometer ORA, sser 12-30

€115.00
plus VAT &
Shipping

Product Images



Description

- The models of the KERN ORA-B series are universal, maintenance-free analog handheld refractometers
- The handy and robust design allows an easy, efficient and durable use in everyday life
- The manual conversion effort is avoided by several selectable scales and excludes application errors
- These scales are specially developed, precisely calculated and verified. They are also characterized by very thin and clear lines
- The optical system and the prism cover are made of special materials, which allow a low-tolerance measurement
- All models are equipped with an eyepiece that can be easily and smoothly adjusted to different visual acuities
- Models marked with "ATC" feature automatic temperature compensation, which enables accurate measurements at different ambient temperatures

Included in delivery:

- Storage box
- Calibration solution
- calibration block, if necessary
- Pipette
- Screwdriver
- Cleaning cloth
- Further accessories are optionally available

Application range honey

The following models are particularly suitable for measuring the Brix value, the water content in honey and the Baumé degree (°Bé) for determining the relative density of liquids

Main application areas:

- Beekeeping
- Honey production

Additional Information

No.	KP-5086
Manufacturer (Brand)	KERN
EAN	4058072146399
Transport temperature	Room temperature
Refractometer measures	Water content
Manufacturer	CORE
Refractometer type	Hand refractometer

