




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN Refractometer ORT,000-1,7001

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- The models of the KERN ORT-1 series are universal, analog Abbe refractometers
- The handy and robust design allows easy, efficient and durable use in everyday life
- The integrated scale allows the use in many applications and offers the best possible safety to read the measurement results accurately

### Included in delivery:

- Calibration solution
- Calibration block
- Pipette
- Screwdriver
- Cleaning cloth
- digital thermometer
- Further accessories are optionally available

### Technical data:

- Overall dimensions W×D×H 180×90×240 mm
- Net weight approx. 1,95 kg

### Field of application: industry/pharmacy/laboratory

The following model is a simple but extremely reliable Abbe refractometer with thermometer. Liquid, solid as well as pasty samples can be evaluated. This refractometer is characterized by its robustness and easy handling. An attractive aluminum transport and storage box is also available as an option. The refractive index  $n_D$  is measured.

### Main application areas:

- Sugar industry (cane sugar)
- Pharmacy
- Beverage industry
- Food industry
- Chemical industry
- Oil industry/refineries
- Laboratories
- Training

## Additional Information

---

No.	KP-5122
Manufacturer (Brand)	KERN
VGKL number	105837801
EAN	4058072146740
Transport temperature	Room temperature
Manufacturer	CORE
Refractometer type	Hand refractometer

