



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Elma ultrasonic unit 12,75 l, 2x600 HF, with heating

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

Robust ultrasonic equipment for efficient cleaning of glass parts, metal parts, cutlery and plastic items.

- Stainless steel oscillating tank
- Ultrasonic power transducer 37 kHz
- optimal sound field distribution due to switchable frequency modulation
- Display of set and remaining time
- LED temperature display
- Timer and continuous operation
- Degas function

## Additional Information

No.	GH-0009
Manufacturer (Brand)	Elma
old neoLab article no.	2-4148
EAN	4058072670078
Transport temperature	Room temperature
Temp MAX	80 °C
Temp MIN	30 °C
Volume	12750 ml
Wide	278 mm
Height	321 mm
Length	365 mm
Depth	365 mm
Operation via	Knobs/Buttons
Weight	7.5 kg
Voltage	240 V
Cleaning equipment and sterilizers Features	with timer with heating with continuous function closure/lid available as accessory
Frequency	37 KHz
Power	1000 W
Material bath	Stainless steel
Bath depth	200 mm
Bath opening	240 x 300 mm

