



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Elma ultrasonic device Mini, 0.8 l, 100 HF

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Robust ultrasonic equipment for efficient cleaning of glass parts, metal parts, cutlery and plastic items.

- Stainless steel oscillating tank
- Ultrasonic power transducer 37 kHz
- optimal sound field distribution due to switchable frequency modulation
- Display of set and remaining time
- LED temperature display
- Timer and continuous operation
- Degas function

Additional Information

No.	GH-0040
Manufacturer (Brand)	Elma
old neoLab article no.	2-4152
EAN	4058072260880
Transport temperature	Room temperature
Volume	800 ml
Wide	116 mm
Height	178 mm
Length	206 mm
Depth	206 mm
Operation via	Knobs/Buttons
Weight	2 kg
Voltage	240 V
Cleaning equipment and sterilizers Features	with timer with continuous function lid/lid available as accessory
Frequency	37 KHz
Power	30 W
Material bath	Stainless steel
Bath depth	60 mm
Bath opening	85 x 190 mm

