



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® ultrasonic unit, 750 ml

€86.80
plus VAT &
Shipping

Product Images



Description

Compact ultrasonic cleaner for cleaning glasses, instruments and other small parts.

- with stainless steel tank and plastic basket
- max. filling capacity: 750 ml
- with timer function
- Frequency: 40 KHz

Additionally with CD holder.

Additional Information

No.	2-3199
Manufacturer (Brand)	neoLab
EAN	4058072090234
Transport temperature	Room temperature
Volume	750 ml
Wide	200 mm
Height	145 mm
Length	175 mm
Depth	175 mm
Cleaning equipment and sterilizers Features	with closure/lid with timer
Bath depth	50 mm
Bath opening	130 x 150 mm
Operation via	Membrane keypad
Frequency	40 KHz
Weight	1.2 kg
Power	40 W
Material bath	Stainless steel
Voltage	230 V

