




neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® Ultrasonic unit, small, 0.5 l

**€236.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Small digital ultrasonic cleaner for safety glasses and other sensitive small items such as tweezers, drills, scissors, fine mechanical parts, micro components, etc..

With removable tank and timer.

Volume: 0.5 l

## Additional Information

---

No.	DV-0738
Manufacturer (Brand)	neoLab
old neoLab article no.	DV-0685
EAN	4066292558465
Transport temperature	Room temperature
Volume	500 ml
Wide	100 mm
Height	124 mm
Length	190 mm
Depth	190 mm
Cleaning equipment and sterilizers Features	with closure/lid with timer
Bath depth	40 mm
Bath opening	80 x 175 mm
Operation via	Membrane keypad
Frequency	42 KHz
Weight	1.4 kg
Power	50 W
Material bath	Stainless steel
Voltage	230 V

