



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Bühler compact shaker KS 15 B, max 15 kg,
shaking surface 400 x 300mm, shaking motion
horizontal**

**€1,997.00
plus VAT &
Shipping**

Product Images



Description

Description Reliable and robust in daily use Scope of delivery Basic unit without vibrating plate, without attachment
Technical data Shaking motion: back and forth (horizontal) Shaking surface: 400 x 300 mm Max. Load: 15 kg Speed control range: 30 - 420 min⁻¹, infinitely variable Stroke: 17 mm Timer: 0 - 120 min / continuous operation Electrical connection: 230 V or 115 V, 50 / 60 Hz (please specify when ordering) Protection class: IP 21 Heat output: approx. 5 - 10 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % External dimensions (W x D x H): 510 x 490 x 150 mm Weight: 18 kg Loading possibilities (examples) Erlenmeyer flasks 100 ml 20 pieces Erlenmeyer flasks 250 ml 12 pieces Erlenmeyer flasks 500 ml 6 pieces Erlenmeyer flasks 1000 ml 4 pieces Separating funnels 250 ml 3 pieces

Additional Information

No.	NB-0143
Manufacturer (Brand)	Bühler
VGKL number	526161730
EAN	4058072625054
Gross weight	22
Transport temperature	Room temperature
Net weight	18.00000 kg
Speed MAX	420 rpm
Country of origin	Germany
Region of origin	Baden-Württemberg
Wide	510 mm
Customs tariff number	84798200
Height	150 mm
Depth	490 mm
Voltage	230 V
for ambient temperature	50 °C
Amplitude	17 mm
Speed MIN	30 rpm
Width tray	400 mm
loadable up to	15 kg
Length shelf	300 mm

