



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Thermo Shaker for 2 microtiter plates

€1,300.00
plus VAT &
Shipping

Product Images



Description

Thermal shaker with circular movements for 2 microtiter plates. Dual heating system for lid and shaker plate to avoid condensation. High precision temperature monitoring and control. Suitable also for very small volume samples. Digital control of temperature, rpm and time. The instrument has a memory for user settings, as well as a special anti-vibration system. It is equipped with a maintenance-free motor for continuous operation.

Additional Information

No.	7-0951
Manufacturer (Brand)	neoLabLine
EAN	4058072091293
Transport temperature	Room temperature
Speed MAX	1300 rpm
Color	yellow
Wide	280 mm
Height	120 mm
Length	305 mm
Depth	305 mm
suitable for	Microtiter plates
Operation via	Membrane keypad
Weight	6.5 kg
Voltage	230 V
for ambient temperature	10 to 40 °C
Shaker properties	with infinitely variable speed control with timer with speed monitoring with continuous function
Amplitude	2 mm
Movement	circling
Speed MIN	100 rpm
Power	140 W
Material block	Aluminum
Temp. accuracy	+/- 0.1 °C
Spare part	7-0250
Temp MIN above ambient temperature	3 °C
Shaker type	Thermal shaker

