



neoLabLine neoMix without cooling function

Product Images



Description

neoMix without cooling function for temperature control and mixing tasks in the laboratory. Exchangeable reaction blocks suitable for vessels from 0.2 to 50 ml or for multiwell plates.

Temperature range in each case up to +100°C, with block temperature homogeneity at < +0.5°C.

The device has a backlit display with plain text display.

Delivery WITHOUT temperature control block !

A temperature control block (7-0922 to 7-0934) is required as an accessory for operation of the thermal shaker.

No.	7-0921
Manufacturer (Brand)	neoLabLine
EAN	4058072162238
Speed MAX	1500 U/Min
Color	gelb
sterile	No
Wide	280 mm
Height	170 mm
Length	220 mm
Depth	220 mm
Operation via	Folientastatur
Weight	8.5 kg
Voltage	230 V
for ambient temperature	5 bis 40 °C
Shaker properties	programmierbar mit stufenloser Drehzahlregelung mit Timer mit Shortmix-Funktion mit Drehzahlüberwachung mit Dauerfunktion
Amplitude	2 mm
Resolution °C	1 °C
Movement	vibrierend kreisend
Speed MIN	200 U/Min
Number of speed levels	stufenlos
Power	85 W
Temp. accuracy	+/- 0.3 °C
Heating rate	9 °C/ Min
Maximum relative humidity at 20 °C	70.0 %
Temp MIN below ambient temperature	5 °C
Shaker type	Thermoschüttler

