



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Heidolph Multi Reax -- EU cable

**€1,595.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The versatile one with 3 mm shaker orbit High sample throughput is one of the everyday challenges faced by analytical laboratories. Often, several analyses have to be prepared at the same time. To simplify these processes and the general workflow, a multi-vortexer has been developed. This can accommodate between 12 ( $\varnothing = 16-32$  mm) and 26 ( $\varnothing = 10-16$  mm) tubes (Falcon tubes or Greiner tubes). Reaction vessels with snap-on lids ("Eppis" or Eppendorff Cups) also find a secure hold. The 3 mm shaker orbit and a speed control of 150-2,000 rpm guarantee excellent mixing results. The integrated timer (up to 999 minutes) allows automatic termination of the shaking function. After the time has expired, a signal tone is emitted. In this way, one can be sure to continue work processes continuously even when completing several tasks at the same time.

## Additional Information

---

No.	BH-1797
Manufacturer (Brand)	Heidolph
VGKL number	524145100
EAN	4058072779450
Gross weight	12.2
Transport temperature	Room temperature
Net weight	10.70000 kg
Speed MAX	2000 rpm
Country of origin	Germany
Region of origin	Bavaria
Packing width	0.38000 m
Packing height	0.29000 m
Packing depth	0.47000 m
Packaging volume	0.05200000 ccm
Customs tariff number	84798200
Speed MIN	150 rpm

