



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph test tube shaker Reax control

€473.00
plus VAT &
Shipping

Product Images



Description

Test tube shaker with strong, circular vibrating motion.

High stability due to heavy die-cast housing.

With switch for interval and continuous operation. In interval operation, the unit switches on when a vessel is pressed onto the receiving plate.

Reax top: with speed control

Reax control: with speed regulation, even at low speeds max. deviation $\pm 2\%$

Technical data:

Shaking frequency: 0-2500/min

Weight: 3.3 kg

Protection class : IP 22

Receiving plate: jars 20 mm

Connection: 230 Volt, 50 Hz, 42 Watt

Additional Information

No.	BH-1725
Manufacturer (Brand)	Heidolph
EAN	4058072195977
Transport temperature	Room temperature
Speed MAX	2500 rpm
Color	gray
Wide	134 mm
Height	105 mm
Length	172 mm
Depth	172 mm
Operation via	Knobs/Buttons
Weight	2.8 kg
Voltage	230 V
for ambient temperature	5 to 40 °C
Shaker properties	maintenance-free with infinitely variable speed control with automatic speed control with press-to-start function with speed monitoring with continuous function including attachment
Figure	2
Amplitude	5 mm
Attachment is suitable for	Reaction vessels. Test tubes. Centrifuge tubes (standard attachment)
Attachment is suitable for type vessel	Reaction vessels cryogenic tubes centrifuge tubes test tubes/tubes
Movement	circling vibrating
Speed MIN	0 rpm
Capacity essay	1
Layer tube	vertical
Power	5 W
Maximum relative humidity at 20 °C	85.0 %

Protection class

IP 22

Shaker type

Vibratory shaker/ vortexer

