




neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® multipurpose vortexer with multifunctional attachment

€380.00
plus VAT &
Shipping

Product Images



Description

With the neoLab multipurpose vortexer, you can mix the standard vessels commonly used in the laboratory with a single attachment.

The universal attachment can be easily removed for cleaning purposes. Three speeds ensure intensive mixing. Press-to-mix and continuous function. Housing made of ABS. Non-slip stand due to rubber feet.

- 1 x attachment Press- to-Mix surface for vortexing
- 1 x attachment for MTP (mounted with rubber bands)
- 2 x holders for 1 x 50 ml tube each (mounted horizontally)
- 2 x holders for 1 x 15 ml tube each (mounted horizontally)
- 2 x holders for 1 x 5 ml tube each (mounted horizontally)
- 2 x holders for 2 x 1.5/2 ml tubes each (mounted vertically)
- 2 x holders for 3 x 0.5 ml tubes each (mounted vertically)
- 2 x holders for 12 x 0.2 ml PCR tubes or strips each (vertically mounted)

Technical data:

Speed: 1000, 2000, 3000 rpm

Orbit: 3.7 mm

Weight: 3.93 kg

Dimensions (LxWxH): 171 x 190 x 197 mm

Additional Information

No.	7-2025
Manufacturer (Brand)	neoLab
EAN	4058072038267
Transport temperature	Room temperature
Speed MAX	3000 rpm
Color	blue
Wide	190 mm
Height	197 mm
Length	171 mm
Depth	171 mm
Operation via	Knobs/Buttons
Weight	3.93 kg
Voltage	240 V
for ambient temperature	2 to 40 °C
Shaker properties	maintenance-free with press-to-start function with continuous function including attachment
Amplitude	3.7 mm
Number of speed levels	3
Attachment is suitable for	PCR tubes. 12 x 0.2 ml (2x) reaction tubes. 1 x 5.0 ml (2x) reaction tubes. 3 x 0.5 ml (2x) reaction vessels. 2 x 1.5/2.0 ml (2x) centrifuge tubes. 1 x 15 ml (2x) centrifuge tubes. 1 x 50 ml (2x) deepwell plate. 96 microtiter plates
Attachment is suitable for type vessel	Reaction vessels cryogenic tubes centrifuge tubes test tubes/tubes
Movement	vibrating
Speed MIN	1000 rpm
Layer tube	horizontal vertical
Power	40 W
Maximum relative humidity at 20 °C	95.0 %
Model Shaker/ Rotator	neoLab multipurpose vortexer (7-2025)

