



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Vortex shaker, continuously adjustable 0-3150 rpm

€270.00
plus VAT &
Shipping

Product Images



Description

The **neolab Vortex Mixer** is a powerful test tube mixer with continuously adjustable speed from 0 to 3,150 rpm. The mixing intensity can be influenced by applying appropriate pressure to the platform. Its sturdy steel housing and non-slip, insulating rubber feet ensure a stable stand and low vibration development during vortexing. The neolab Vortex Mixer can be operated with press-to-mix or continuous function.

The range of applications can also be usefully extended by a large selection of accessory attachments.

Additional Information

No.	7-2020
Manufacturer (Brand)	neoLab
EAN	4058072038236
eCl@ss Article Code	32012608
Transport temperature	Room temperature
Speed MAX	3200 rpm
Color	yellow
Wide	130 mm
Height	180 mm
Length	165 mm
Depth	165 mm
Operation via	Knobs/Buttons
Weight	5 kg
Voltage	230 V
Shaker properties	maintenance-free with infinitely variable speed control with press-to-start function with continuous function including attachment (Further) attachments are available as accessories
Amplitude	2.8 mm
Movement	vibrating
Speed MIN	0 rpm
Number of speed levels	stepless
Power	20 W
Attachment is suitable for	Reaction vessels. Test tubes. Centrifuge tubes (standard attachment)
Spare part	7-2021 7-2022
Model Shaker/ Rotator	neoLab test tube mixer Vortex (7-2020)
Shaker type	Vibratory shaker/ vortexer

