



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Vortex Genie 2, 2700 rpm

**€495.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Vortex-Genie®, the universal small shaker with extra power!

High-speed vortexer.

Fast vortexing of problematic samples, easy to use. Circular motion, 6 mm orbit. Each including rubber plate (inside 2.5 cm Ø) and with 8.6 cm platform. Attachments easily exchangeable.

Order No. 7-0092 VortexGenie®2

With press-to-mix function

Weight: 4 kg

Power supply: 240 V, 50 Hz

## Additional Information

No.	7-0092
Manufacturer (Brand)	neoLab
EAN	4058072027018
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	2700 rpm
Color	blue
sterile	No
Wide	122 mm
Height	165 mm
Length	165 mm
Depth	165 mm
Operation via	Knobs/Buttons
Weight	4 kg
Voltage	230 V
for ambient temperature	0 to 40 °C
Shaker properties	with infinitely variable speed control with press-to-start function with continuous function including attachment (Further) attachments are available as accessories
Amplitude	1 mm
Figure	1
Movement	vibrating circulating
Speed MIN	600 rpm
Number of speed levels	stepless
Power	80 W
Attachment is suitable for	Reaction vessels. Test tubes. Centrifuge tubes (standard attachment)
Attachment is suitable for type vessel	Reaction vessels cryogenic tubes centrifuge tubes test tubes/tubes

Spare part	7-0510 7-0513 7-0514
Maximum relative humidity at 20 °C	95.0 %
Model Shaker/ Rotator	neoLab Vortex Mixer Genie 2 (7-0092)
Shaker type	Vibratory shaker/ vortexer

